

POLYMED[®]
by **Sempermed**[®] *STERILE*
Latex examination gloves

Industrial

Exam

Surgical

Latex • Sterile • Powder-Free



Sempermed[®]
We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



The characteristics that make Polymed® a leading choice in hand protection are now available in sterile pairs. Polymed® Sterile is designed for critical environments that require a sterile glove, but not necessarily a more expensive surgical glove. Polymed® Sterile can be especially effective for wound care, trauma, and burn care. Sterile exam gloves are also appropriate for OB/GYN, CCU, and ER departments. The gloves are powder-free and low protein.

FEATURES

- 50 micrograms or less of total water-extractable protein per gram*
- Powder-Free
- Sterile

BENEFITS

- Reduced potential for allergic reaction
- Eliminates the potential aerosolization of glove powder
- For aseptic environments

GLOVE COUNT

40 Pair/Box

6 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.08 mm/3.2 mils min.

Palm: 0.08 mm/3.2 mils min.

PHYSICAL PROPERTIES*

Before aging:

Tensile Strength: 18 MPa min.

Elongation: 650% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

PRODUCT COMPOSITION

- Natural Rubber Latex
- Sulfur, Zinc Oxide, Organic Accelerators (dithiocarbamates, benzothiazolate)

STERILIZATION

- Gamma Radiation
- Minimum Dose: 25 kGy
- Indicator Label on Dispenser Box and Carton

Caution: This product contains natural rubber latex which may cause allergic reactions. Safe use of this product by or on latex sensitized individuals has not been established.

Order #	Size	Length (min.)	Palm Width (min.)
PSE102	S	220mm	80mm +/- 10mm
PSE103	M	230mm	95mm +/- 10mm
PSE104	L	230mm	111mm +/- 10mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558