

SemperCARE[®] CR

by Sempermed[®] *STERILE*

Chloroprene examination gloves

Industrial

Exam

Surgical

Non-Latex • Sterile • Powder-Free



Sempermed[®]
We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



SemperCare® CR Sterile is designed for critical environments that require a sterile glove, but not necessarily a more expensive surgical glove. SemperCare® CR Sterile singles are especially appropriate for Labor & Delivery and OB/GYN.

The gloves are 100% latex-free and are manufactured using an advanced polychloroprene formula. These powder-free gloves offer the fit and comfort you expect from the finest latex.

FEATURES

- No natural rubber latex
- Powder-Free
- Look and feel of natural rubber latex
- Sterile

BENEFITS

- A latex alternative
- Eliminates the potential aerosolization of glove powder
- Polychloroprene formulation provides fit and comfort of natural rubber latex, without protein allergens
- For aseptic environments

GLOVE COUNT

80 Singles/Box

6 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.05 mm/2.0 mils min.

Palm: 0.05 mm/2.0 mils min.

PHYSICAL PROPERTIES*

Before Aging:

Tensile Strength: 14 MPa min.

Elongation: 500% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

PRODUCT COMPOSITION

- Synthetic Rubber Elastomers
- Zinc Oxide
- Organic accelerator
- Hindered Polymeric Phenol
- Sulfur
- Titanium Dioxide/Bentonite
- Potassium Hydroxide

STERILIZATION

- Gamma Irradiation min. 25 kGy
- Indicator Label on Dispenser Box and Carton

Order #	Size	Length (min.)	Palm Width (avg.)
CRS103	M	230mm +/- 5mm	95mm +/- 10mm
CRS104	L	230mm +/- 5mm	110mm +/- 10mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558

