

# ***SemperGuard<sup>®</sup> Latex***

**INDUSTRIAL GLOVES**

Surgical

Industrial

Exam

Latex • Powder-Free and Powdered



**Sempermed<sup>®</sup>**

We are people protection.™

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



SemperGuard® Latex gloves are designed to provide the right combination of protection, sensitivity, economy, and value. Our latex gloves are strong and elastic providing maximum protection and comfort. These gloves are available in powder-free and powdered.

## FEATURES

- Flexible and elastic
- Enhanced puncture resistance
- Meets FDA CFR 21 part 177 requirements for food handling

## BENEFITS

- Enhanced comfort and fit
- Reduced possibility of tears and holes
- Safe for food contact applications

### GLOVE COUNT (By Weight)

100 Gloves/Box (90 Gloves/Box XL)  
10 Boxes/Case

### THICKNESS\* (Single Wall)

Finger: 0.08 mm/3.2 mils min.  
Palm: 0.08 mm/3.2 mils min.

### PHYSICAL PROPERTIES\*

Before aging:

Tensile Strength: 18 MPa min.  
Elongation: 650% min.

### SURFACE TREATMENT

Powdered with U.S.P Absorbable Dusting Powder

### PRODUCT CHARACTERISTICS

- Natural Rubber Latex
- Ambidextrous
- Beaded Cuff
- Single Use Only
- Absorbable Dusting Powder U.S.P. (SemperGuard® Latex Powdered)

**CAUTION:** This product contains natural rubber latex which may cause allergic reactions.

\*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request. All Sempermed gloves are single use only.

### POWDER-FREE

Order #	Size	Length (min.)	Palm Width (min.)
INDPFT102	S	230mm	85mm +/- 5mm
INDPFT103	M	230mm	95mm +/- 5mm
INDPFT104	L	230mm	112mm +/- 5mm
INDPFT105	XL	230mm	115mm +/- 5mm

### POWDERED

Order #	Size	Length (min.)	Palm Width (min.)
INDPS102	S	220mm	83mm +/- 5mm
INDPS103	M	230mm	95mm +/- 5mm
INDPS104	L	230mm	105mm +/- 5mm
INDPS105	XL	230mm	117mm +/- 5mm

**Sempermed®**

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762  
Phone: 1.800.366.9545 • Fax: 1.800.763.5491  
Outside the U.S.: 727.787.7250 • Fax: 727.787.7558