

# SEMPERSHIELD™

EXTENDED CUFF

## NITRILE EXAMINATION GLOVES

Industrial

Exam

Surgical

Non-Latex • Powder-Free • Textured



**Sempermed®**   
We are people protection.®

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



With extra thickness and an extended cuff, SemperShield™ Nitrile powder-free exam gloves are designed for high-risk applications where your safety and protection is a must. A textured surface promises a secure grip without compromising your sensitivity. SemperShield™ Nitrile offers all the benefits of a nitrile glove plus it's safe for use with chemotherapy drugs.

## FEATURES

- Extra thickness and an extended cuff
- Excellent chemical resistance
- Tested and approved for use with chemotherapy drugs
- Manufactured without natural rubber
- Powder-Free

## BENEFITS

- Additional barrier protection for the wrist and lower forearm
- Protects you against a variety of chemicals
- Safe for use with chemotherapy drugs
- Eliminates exposure to latex
- Eliminates aerosolization of glove powder

### GLOVE COUNT (By Weight)

50 Gloves/Box  
10 Boxes/Case

### THICKNESS\* (Single Wall)

Finger: 0.05 mm/2.0 mils min.  
Palm: 0.05 mm/2.0 mils min.

### PHYSICAL PROPERTIES\*

Before aging:

Tensile Strength: 14 MPa min.  
Elongation: 500% min.



Manufactured in accordance with the National Fire Protection Association standard on protective clothing for emergency medical operations.

\*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.  
All Sempermed gloves are single use only.

### PRODUCT CHARACTERISTICS

- Synthetic Nitrile Rubber (NBR)
- Textured
- Powder-Free
- Ambidextrous
- Beaded Cuff
- Single Use Only

Chemotherapy drugs tested: Carmustine, Cisplatin, Cyclophosphamide, Dacarbazine, Doxorubicin Hydrochloride, Etoposide, 5-Fluorouracil, Paclitaxel, Thio-Tepa (call 800.366.9545 for copy of test data).

Order #	Size	Length (min.)	Palm Width (avg.)
SSNF102	S	290mm	80mm +/- 10mm
SSNF103	M	290mm	95mm +/- 10mm
SSNF104	L	290mm	110mm +/- 10mm
SSNF105	XL	290mm	120mm +/- 10mm
SSNF106	XXL	290mm	122mm +/- 5mm



### HOW WE GLOVE THE PLANET WITH SEMPERSHIELD® NITRILE

- ✓ Inner boxes made with 80% recycled material
- ✓ Outer cartons made with 90% recycled material

**Sempermed**®

We are people protection.®

www.SempermedUSA.com

13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558



Glove the Planet™ is Sempermed USA Inc.'s commitment to reducing our environmental impact. By thinking ahead and consuming less we're able to do more to minimize our global footprint with every step we take. Our ethos of sustainability drives innovative production and distribution practices that lead to greater efficiencies, resulting in less waste, reduced emissions, and a healthier planet for future generations. To learn more, visit [www.GlovethePlanet.com](http://www.GlovethePlanet.com).



TM/® Trademark/registered trademark of Sempermed USA, Inc. 12/09