

sempermed® 
DERMA plus

Exam

Surgical

Industrial

Latex • Powdered • Micro-Rough



Sempermed® 
We are people protection.®

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



Sempermed® Derma Plus powdered surgical gloves were designed with the surgeon in mind. These gloves deliver comfort, barrier protection, and value. They offer ease of donning for both wet and dry hands. The broader, anatomical design allows for an excellent fit, while the micro-rough surface provides for sensitivity and an effective grip, wet or dry. Beaded cuffs increase strength and protection.

FEATURES

- Broad, anatomical design
- Micro-rough surface
- Beaded cuff
- Size stamped on cuff
- Sterility indicator label on dispenser

BENEFITS

- Exceptional fit
- Provides an effective grip, wet or dry, while maintaining tactile sensitivity
- Increased strength and protection
- Easy to find proper glove size
- Confirmation of sterility

GLOVE COUNT (By Weight)

50 Pair/Box

6 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.10 mm/4.0 mils min.

Palm: 0.10 mm/4.0 mils min.

PHYSICAL PROPERTIES*

Before aging:

Tensile Strength: 24 MPa min.

Elongation: 750% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

SURFACE TREATMENT

Powdered with U.S.P Absorbable Dusting Powder

SURFACE TREATMENT

- Natural Rubber Latex
- Sulfur
- Zinc Oxide
- Organic accelerator (dithiocarbamate)
- Polymeric sterically hindered phenol

CAUTION: This product contains natural rubber latex which may cause allergic reactions.

Order #	Size	Length (min.)	Palm Width (min.)
DP600	6	265mm	76mm +/- 6mm
DP650	6.5	265mm	83mm +/- 6mm
DP700	7	265mm	89mm +/- 6mm
DP750	7.5	265mm	95mm +/- 6mm
DP800	8	265mm	102mm +/- 6mm
DP850	8.5	265mm	108mm +/- 6mm

Sempermed®

We are people protection.®

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762

Phone: 1.800.366.9545 • Fax: 1.800.763.5491

Outside the U.S.: 727.787.7250 • Fax: 727.787.7558