

sempermed[®] 
Surgical gloves **supreme**

Exam

Surgical

Industrial

Latex • Powder-Free • Micro-Rough



Sempermed[®] 
We are people protection.[™]

Sempermed is a vertically integrated manufacturer that provides focused sales, marketing, and logistical support to those in need of hand protection.



Surgeons' hands are the most delicate instruments in the surgical suite. They need to be prepared. And protected. That's why we developed Sempermed® Supreme surgical gloves. These unique gloves are anatomically designed to reduce hand fatigue and increase comfort. They are powder-free to eliminate aerosolization and skin irritation. A special multi-layer copolymer composition provides optimal protection during surgery.

FEATURES

- Powder-Free
- Multi-layer coating
- Broad, anatomical design
- Micro-rough surface
- Beaded cuff
- Size stamped on cuff
- Sterility indicator label on dispenser

BENEFITS

- Eliminated aerosolization of glove powder
- Copolymer design offers comfort and ease of donning, wet or dry
- Exceptional fit
- Provides an effective grip, wet or dry, while maintaining tactile sensitivity
- Increased strength and protection
- Easy to find proper glove size
- Confirmation of sterility

GLOVE COUNT

50 Pair/Box

6 Boxes/Case

THICKNESS* (Single Wall)

Finger: 0.10 mm/4.0 mils min.

Palm: 0.10 mm/4.0 mils min.

PHYSICAL PROPERTIES*

Before aging:

Tensile Strength: 18 MPa min.

Elongation: 750% min.

*The specifications provided reflect the minimum acceptable requirements of ASTM and FDA where applicable. Gloves offered by Sempermed USA, Inc. meet or exceed these physical requirements. Independent laboratory test results are available upon request.

All Sempermed gloves are single use only.

PRODUCT COMPOSITION

- Natural Rubber Latex
- Sulfur
- Zinc Oxide
- Organic Accelerators
- Antioxidants
- Synthetic Polymer

CAUTION: This product contains natural rubber latex which may cause allergic reactions.

Order #	Size	Length (min.)	Palm Width (min.)
SPFP550	5.5	245mm	70mm +/- 6mm
SPFP600	6	265mm	76mm +/- 6mm
SPFP650	6.5	265mm	83mm +/- 6mm
SPFP700	7	265mm	89mm +/- 6mm
SPFP750	7.5	265mm	95mm +/- 6mm
SPFP800	8	265mm	102mm +/- 6mm
SPFP850	8.5	265mm	108mm +/- 6mm
SPFP900	9	265mm	114mm +/- 6mm

Sempermed®

We are people protection.™

www.SempermedUSA.com • 13900 49th Street North • Clearwater, FL 33762
 Phone: 1.800.366.9545 • Fax: 1.800.763.5491
 Outside the U.S.: 727.787.7250 • Fax: 727.787.7558

